

AERODROME GROUND
MOVEMENT CHART - ICAO

APRON ELEV
21 M

SEPANG /
KL INTERNATIONAL
AIRPORT
TERMINAL 1

TABLE 3

| TAXIWAY WIDTH AND BEARING STRENGTH | | | |
|------------------------------------|-------------|-------|------------------------------------|
| NO | APN TAXIWAY | WIDTH | BEARING STRENGTH |
| 62 | G | 29 | PCR 448/F/C/W/U |
| 63 | H | 25 | PCR 448/F/C/W/U |
| | TAXIWAY | WIDTH | BEARING STRENGTH |
| 64 | H3 | 25 | PCR 448/F/C/W/U |
| 65 | H5 | 25 | PCR 448/F/C/W/U |
| | APN TAXIWAY | WIDTH | BEARING STRENGTH |
| 66 | K | 25 | PCR 540/R/C/W/U PCR 448/F/C/W/U |
| 67 | L | 24 | PCR 540/R/C/W/U PCR 448/F/C/W/U |
| 68 | M | 25 | PCR 540/R/C/W/U |
| | TAXIWAY | WIDTH | BEARING STRENGTH |
| 69 | N1 | 50 | PCR 448/F/C/W/U |
| 70 | N2 | 35 | PCR 448/F/C/W/U |
| | APN TAXIWAY | WIDTH | BEARING STRENGTH |
| 71 | S1 | 24.8 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 72 | S2 | 34 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 73 | S3 | 23 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 74 | S4 | 35.8 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 75 | S5 | 23 | PCR 448/F/C/W/U |
| 76 | S6 | 24 | PCR 448/F/C/W/U |
| 77 | S7 | 23 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 78 | S8 | 44.5 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 79 | S9 | 26 | PCR 540/R/C/W/U |
| 80 | S10 | 42.5 | PCR 540/R/C/W/U |
| 81 | T2 | 23 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 82 | T3 | 29 | PCR 448/F/C/W/U |
| 83 | T4 | 32 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 84 | T5 | 30 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 85 | T6 | 24 | PCR 448/F/C/W/U |
| 86 | T7 | 23 | PCR 448/F/C/W/U PCR 540/R/C/W/U |
| 87 | T13 | 26 | PCR 540/R/C/W/U |

CHANGES : TABLE FOR BEARING STRENGTH